

Reax top

P/N: 541-10000-00



The classic with 5 mm shaking orbit

This model is particularly suitable for simple mixing tasks. In short-term operation mode, the shaking motion is triggered by pressure on the test tube tray. This mode is suitable for mixing up solutions, for example.

Persistent cell pellets can be dissolved gently but powerfully in continuous operation. With a speed of up to 3,200 rpm and a shaking orbit of 5 mm, reliable and uniform distribution of solids, and homogenization of highly viscous media are guaranteed.

Reax top - Technical Data

Overheat protection	self-resetting
Permissible ambient conditions	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Weight	2.8 kg
Protection class IEC 60529	IP 22
Timer	-
Rotation speed range	200 - 3,300 rpm
Power input	51 W
Dimensions (w/d/h)	134 x 172 x 105 mm
Max. Load	-
Operating Mode	automatic or continuous
Rotation speed setting	analog
Orbit / Stroke	5 mm
Platform size	-
Accessories included	-
Angle	-
Motion	vibrating